

Monovisc®

Hyaluronic Acid Single Injection

Medical Information Brochure

Structured product education for osteoarthritis pain management, designed for the Platinum Pharma website brochure download section.



Single injection



Up to 6 months*



Ultra-pure HA



HA helps lubricate, cushion,
and mechanically support the joint.



Product overview: Monovisc® is a lightly cross-linked, high molecular weight hyaluronic acid viscosupplement administered as a single intra-articular injection to help relieve symptoms of osteoarthritis by improving the lubricating properties of synovial fluid. Claims and indication details should be verified against the Egypt-approved IFU and local registration documents before publication.

*Pain relief duration and clinical claims must be verified against the Egypt-approved IFU and local registration documents before publication.

Product facts and administration

Key product facts

- Category: hyaluronic acid viscosupplement
- Composition: lightly cross-linked sodium hyaluronate
- Presentation: sterile viscoelastic preparation in a single-use syringe
- Strength: 22 mg/mL in a 4.0 mL syringe
- Source: ultra-pure, non-avian hyaluronic acid produced by bacterial fermentation
- Intended use: osteoarthritis symptom / pain relief according to local approval

Administration summary

- Single intra-articular injection into the selected joint space
- Administered by a qualified healthcare professional using strict aseptic technique
- Remove joint effusion, if present, before injection
- Product selection and suitability should follow the approved IFU and physician judgment
- Post-procedure monitoring and patient counseling should follow local practice

Treatment course visual




Post-injection patient advice

As with any invasive joint procedure, patients are advised to avoid strenuous activity or prolonged weight-bearing activity, such as running or tennis, within 48 hours after the intra-articular injection. Patients should follow their treating physician's instructions.

Understanding osteoarthritis and viscosupplementation

A healthy joint contains synovial fluid, a viscous liquid that lubricates and absorbs shocks. In osteoarthritis, changes in synovial fluid and degeneration of joint structures may contribute to pain, stiffness, reduced mobility, and ongoing discomfort. Hyaluronic acid supplementation is intended to improve the lubricating properties of the joint.



1 Lubrication

Hyaluronic acid supplementation is intended to improve synovial-fluid lubrication and help reduce friction within the joint.

2 Cushioning

The gel-like viscoelastic nature of Monovisc® helps support cushioning during movement and weight-bearing activities.

3 Single-injection convenience

Monovisc® is designed as a single-injection viscosupplement, providing a convenient treatment option within the approved labeling.

What is Monovisc®?

Monovisc® is a sterile, viscoelastic preparation of lightly cross-linked sodium hyaluronate supplied in a single-use syringe. It is used as a hyaluronic acid viscosupplement for osteoarthritis-related joint pain according to the applicable regional approval.



Recommended website brief

Monovisc® is a lightly cross-linked, high molecular weight hyaluronic acid viscosupplement designed to help relieve symptoms of osteoarthritis by supplementing the synovial fluid. It is administered by a qualified healthcare professional as a single intra-articular injection, according to approved product labeling and local registration.

Evidence highlights and safety

Clinical evidence summary

Monovisc® website and brochure materials describe long-lasting pain relief from a single intra-articular injection.

The product website states that a single injection can provide osteoarthritis pain relief lasting up to 6 months, with relief in the knee felt as early as 2 weeks.



*Duration of effect is individual and subject to clinical response, local labeling, and medical approval.



Long-term relief

Website materials describe pain relief lasting up to 6 months, subject to local labeling and medical approval.



Single-injection convenience

A one-injection treatment course may support convenience for appropriate patients and clinicians.



Well-tolerated HA

Monovisc® is presented as an ultra-pure, non-avian hyaluronic acid viscosupplement.

Essential safety and suitability notes

Do not administer to patients with known hypersensitivity to hyaluronate preparations or gram-positive bacterial proteins, or in joints with infection or skin disease at the injection site. Transient pain, swelling, or injection-site reactions may occur after intra-articular injection. Monovisc® should be injected only by medical professionals trained in accepted intra-articular injection technique. Final publication should be checked against the Egypt-approved IFU and local registration file.

Core sources used: Monovisc.ca product page and resources; Anika Monovisc brochure / IFU references; KneePainRelief.ca comparison section. Regulatory note: final publication should be checked against the Egypt-approved IFU/registration file and Platinum Pharma/Anika medical approval.